



AKHYAR

Professional Diagnostic
Instruments



Table of Contents



McIntosh Laryngoscope	Click Here
McIntosh Laryngoscope Fiber Optic	Click Here
Miller Laryngoscope	Click Here
Miller Laryngoscope Fiber Optic	Click Here
McIntosh Integrated Laryngoscope	Click Here
McCoy Flexi Tip Fiber Optic	Click Here
McCoy Flexi Tip Integrated Fiber Optic	Click Here
Mini Otoscope Fiber Optic	Click Here
Ent Devices (Sets)	Click Here
Diagnostic (Sets)	Click Here
Standard Accessories	Click Here
Stethoscopes	Click Here
Different Types of Packing Boxes and Pouches	Click Here



McIntosh Laryngoscope

Mc – Intosh Laryngoscopes Standard Illumination Adult Size



AK-100-05

Handle Standard Illumination
Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
■ Set of 3 Blades Size 0, 1, 2 and Adult Handle	■ AK-200-30
■ Set of 3 Blades Size 1, 2, 3 and Adult Handle	■ AK-200-31
■ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	■ AK-200-40
■ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	■ AK-200-41

Sets available in hard case or soft pouch



AK-100-00

Mc - Intosh Blade Standard Illumination
Size: 0



AK-100-01

Mc - Intosh Blade Standard Illumination
Size: 1



AK-100-02

Mc - Intosh Blade Standard Illumination
Size: 2



AK-100-03

Mc - Intosh Blade Standard Illumination
Size: 3



AK-100-04

Mc - Intosh Blade Standard Illumination
Size: 4

Specification:

- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade

McIntosh Laryngoscope - Fiber Optic



Mc - Intosh Laryngoscopes Fiber Optic Illumination Adult Size



AK-101-09

Handle Fiber Optic Illumination
Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
■ Set of 3 Blades Size 0, 1, 2 and Adult Handle	■ 21-200-30
■ Set of 3 Blades Size 1, 2, 3 and Adult Handle	■ 21-200-31
■ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	■ 21-200-40
■ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	■ 21-200-41

Sets available in hard case or soft pouch



AK-101-04

Mc - Intosh Blade Fiber Optic Illumination
Size: 0



AK-101-05

Mc - Intosh Blade Fiber Optic Illumination
Size: 1



AK-101-06

Mc - Intosh Blade Fiber Optic Illumination
Size: 2



AK-101-07

Mc - Intosh Blade Fiber Optic Illumination
Size: 3




AK-101-08

Mc - Intosh Blade Fiber Optic Illumination
Size: 4

Specification:


- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



Miller Laryngoscope

Miller Laryngoscopes Standard Illumination Child / Infant Size



AK-102-09

Handle Standard Illumination
Size: Child / Infant

Handle options: Adult & Stubby

Standard Sets

Contents	Code
■ Set of 3 Blades Size 00, 0, 1 and Child Handle	■ 22-200-30
■ Set of 3 Blades Size 0, 1, 2 and Child Handle	■ 22-200-31
■ Set of 4 Blades Size 00, 0, 1, 2 and Child Handle	■ 22-200-40
■ Set of 4 Blades Size 0, 1, 2, 3 and Child Handle	■ 22-200-41

Sets available in hard case or soft pouch



AK-102-04

Miller Blade Standard Illumination
Size: 00



AK-102-05

Miller Blade Standard Illumination
Size: 0



AK-102-06

Miller Blade Standard Illumination
Size: 1



AK-102-07

Miller Blade Standard Illumination
Size: 2



AK-102-08

Miller Blade Standard Illumination
Size: 3

Specification:

- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x AA size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade

Miller Laryngoscope - Fiber Optic



Miller Laryngoscopes Fiber Optic Illumination

Adult Size



AK-103-04

Handle Fiber Optic Illumination
Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
■ Set of 3 Blades Size 0, 1, 2 and Adult Handle	■ 23-200-30
■ Set of 3 Blades Size 1, 2, 3 and Adult Handle	■ 23-200-31
■ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	■ 23-200-40
■ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	■ 23-200-41

Sets available in hard case or soft pouch



AK-103-04

Miller Blade Fiber Optic Illumination
Size: 0



AK-103-05

Miller Blade Fiber Optic Illumination
Size: 1



AK-103-06

Miller Blade Fiber Optic Illumination
Size: 2



AK-103-07

Miller Blade Fiber Optic Illumination
Size: 3



AK-103-08

Miller Blade Fiber Optic Illumination
Size: 4

Specification:

- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



McInsosh Integrated Laryngoscope

Mc – Intosh Integrated Laryngoscopes Fiber Optic Illumination Adult Size



AK-104-10

Handle Fiber Optic Illumination
Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
■ Set of 3 Blades Size 00, 0, 1 and Adult Handle	■ 24-200-30
■ Set of 3 Blades Size 1, 2, 3 and Adult Handle	■ 24-200-31
■ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	■ 24-200-40
■ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	■ 24-200-41

Sets available in hard case or soft pouch



AK-104-04

Mc - Intosh Blade Fiber Optic Illumination
Size: 00



AK-104-05

Mc - Intosh Blade Fiber Optic Illumination
Size: 0



AK-104-06

Mc - Intosh Blade Fiber Optic Illumination
Size: 1



AK-104-07

Mc - Intosh Blade Fiber Optic Illumination
Size: 2



AK-104-08

Mc - Intosh Blade Fiber Optic Illumination
Size: 3




AK-104-09

Mc - Intosh Blade Fiber Optic Illumination
Size: 4

Specification:


- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade

McCoy Flexi Tip - Fiber Optic



McCoy Flexi - Tip Laryngoscopes Fiber Optic Illumination Adult Size



AK-105-04

Handle Fiber Optic Illumination
Size: Adult



AK-105-04

McCoy Blade Fiber Optic Illumination
Size: 1



AK-105-04

McCoy Blade Fiber Optic Illumination
Size: 2



AK-105-04

McCoy Blade Fiber Optic Illumination
Size: 3

Standard Sets


Contents	Code
■ Set of 2 Blades Size 2, 3 and Adult Handle	■ 25-200-30
■ Set of 2 Blades Size 3, 4 and Adult Handle	■ 25-200-31
■ Set of 3 Blades Size 2, 3, 4 and Adult Handle	■ 25-200-40

Sets available in hard case or soft pouch



Specification:

- Flexible Tip for difficult Intubation
- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



McCoy Flexi Tip Integrated

McCoy Flexi - Tip Integrated Laryngoscopes Fiber Optic Illumination Adult Size



AK-106-07

Handle Fiber Optic Illumination
Size: Adult



AK-106-04

McCoy Blade Fiber Optic Illumination
Size: 1



AK-106-05

McCoy Blade Fiber Optic Illumination
Size: 2



AK-106-06

McCoy Blade Fiber Optic Illumination
Size: 3

Standard Sets

Contents	Code
Set of 2 Blades Size 2, 3 and Adult Handle	26-200-30
Set of 2 Blades Size 3, 4 and Adult Handle	26-200-31
Set of 3 Blades Size 2, 3, 4 and Adult Handle	26-200-40

Sets available in hard case or soft pouch



Specification:

- Flexible Tip for difficult Intubation
- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade

Mini Otoscope Fiber Optic



AK-300-04

Mini Otoscope Fiber Optic
Complete Set of Disposable Cannulas



A

B

C

D

Color Options

Specification:

- Manufactured From Plastic Material
- 3x Magnifying glass
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x AA size batteries
- Soft Pouch available in matching colors
- Fiber Optic Light Source



Ent Devices - Fiber Optic

ENT Devices Fiber Optic Illumination



AK-200-04
Brass Otoscope Fiber Optic
Brass Otoscope Conventional



AK-200-05
ENT Set Complete



AK-200-06
Brass Ophthalmoscope

Specification:

- Manufactured From Brass Chrome plated
- 3x Magnifying glass
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available
- Fiber Optic Light Source

Different Diagnostic Sets



Fiber Optic Diagnostic Set
AK-2000-02



Small Diagnostic Set
AK-2000-03



Universal Diagnostic Set
AK-2000-04



Standard Accessories

Standard Accessories

Nasal Speculum



AK-400-04



AK-400-05



AK-400-06

Bulbs



AK-500-05



AK-500-06



AK-500-07



AK-500-08

Ear Speculas



AK-600-05



AK-600-06



AK-600-07



AK-600-08

Ear Speculas & Testing Tube



AK-600-15



AK-600-20



AK-600-30



Baby Prestige stethoscope
AK-1000-01



Child Prestige stethoscope
AK-1000-02



Kirchner Colorscope stethoscope
AK-1000-03



Standard stethoscope
AK-1000-04



Dual Head stethoscope
AK-1000-05



Different Types of Packing Boxes

Packaging Options



AK-800-05



AK-800-06



AK-800-07



AK-600-08



Akhyar Surgical Professionals

Kashmir Road Street Aslam Nagar
Pacca Graha Town Sialkot 51310
Pakistan

Phone: +92 332 04444 57

E-Mail: info@aksurgicalpro.com

Website: www.aksurgicalpro.com