


Liposuction Cannula Catalogue

Manufacturers & Exporters



AKHYAR

Liposuction Cannula



www.aksurgicalpro.com



BLACK HANDLE

Product Code:AK-LP-10

Product Description:

Black Handle

Handle made of aluminum with type 2 black anodized, 6.65" length 1" diameter



BLUE HANDLE

Product Code:AK-LP-15

Product Description:

BLUE Handle

Handle made of aluminum with type 2 blue anodized, 6.65" length 1" diameter



GOLD HANDLE

Product Code:AK-LP-20

Product Description:

Gold Handle

Handle made of aluminum with type 2 gold anodized, 6.65" length 1" diameter



HANDLE w/SWIVEL BARB

Product Code:AK-LP-25

Product Code:

Product Description:

Handle made of aluminum with type 2 blue anodized, 6.65" length, 1" diameter
(for Large Handle, 3/4" for Small Handle)



HANDLE

Product Code:AK-LP-30

Product Code:AK-LP-35

Product Description:

Handle made of aluminum with type 2 blue anodized, 6.65" length 1" diameter
(for Large Handle, 3/4" for Small Handle)



Liposuction Cannula Catalogue

Liposuction Cannula Handles Cont.



STAINLESS STEEL TOOMEY 60cc HUB

Product Code:AK-LP-50

Product Description:

Stainless Steel and Stainless Steel hub



ALUMINUM TOOMEY 60cc HUB

Product Code:AK-LP-55

Product Description:

Hub made of Aluminum with type 2 blue anodized, 1.55" length 1.31" diameter



QUICK SCREW HUB

Product Code:AK-LP-60

Product Description:

Hub made of Stainless Steel, 2.1" length 7/8" diameter



SMALL QUICK DISCONNECT HUB

Product Code:AK-LP-65

Product Description:

Hub made of Stainless Steel, 1.93" length 5/8" diameter



QUICK DISCONNECT HUB

Product Code:AK-LP-70

Product Description:

Hub made of Stainless Steel, 1.6" length 7/8" diameter



QUICK LOCK HUB

Product Code:AK-LP-75

Product Description:

Hub made of Stainless Steel, 1.65" length 1/2" diameter



LUER LOCK HUB

Product Code:AK-LP-80

Product Description:

Luer Lock Hub

Hub made of Stainless Steel, 1.16" length 3/8" diameter



MERCEDES

Product Code:AK-LP-80

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



TRI PORT

Product Code:AK-LP-85

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



LAS VEGAS

Product Code:AK-LP-90

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



STANDARD 3 PORT

Product Code:AK-LP-95

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



STANDARD 2 PORT

Product Code:AK-LP-100

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



STANDARD SINGLE PORT

Product Code:AK-LP-105

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



3 PORT 60 DEGREE INDEX

Product Code:AK-LP-110

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



3 PORT 60 DEGREE INDEX

Product Code:AK-LP-115

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



BASKET

Product Code:AK-LP-120

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



SPATULA

Product Code:AK-LP-125

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



Liposuction Cannula Catalogue

Cannula Hole Pattern & Cont.



SPATULA 2 PORT

Product Code:AK-LP-130

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



SPATULA 3 PORT

Product Code:AK-LP-135

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



SPATULA 3 PORT

Product Code:AK-LP-140

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



COBRA

Product Code:AK-LP-145

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



FUENTES

Product Code:AK-LP-150

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



SAYLAN

Product Code:AK-LP-155

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



CANDY CANE

Product Code:AK-LP-160

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



MODERN SPIRAL

Product Code:AK-LP-165

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



TOLEDO SINGLE PORT

Product Code:AK-LP-170

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



GARDEN SPRAY ANESTHESIA NEEDLE

Product Code:AK-LP-175

Available Diameter & Length

Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



MULTI HOLE SPRAY NEEDLE

Product Code:AK-LP-180

Available Diameter & Length

Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



BODY ANESTHESIA NEEDLE

Product Code:AK-LP-185

Available Diameter & Length

Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



Infusion 6 hole Needle

Product Code:AK-LP-190

Available Diameter & Length

Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



Infusion 8 hole Needle

Product Code:AK-LP-195

Available Diameter & Length

Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



MICRO INJECTION NEEDLE, STYLE 3

Product Code:AK-LP-200

Available Diameter & Length

Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm.



MICRO INJECTION NEEDLE, STYLE 2

Product Code:AK-LP-205

Available Diameter & Length

Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm.



MICRO INJECTION NEEDLE, STYLE 1

Product Code:AK-LP-210

Available Diameter & Length

Diameter: 18ga,17ga,16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



Liposuction Cannula Catalogue

Reinjection Micro Cannula & Cont.



MICRO INJECTION NEEDLE "V" DISSECTOR

Product Code:AK-LP-215

Available Diameter & Length
Diameter:18ga,17ga,16ga.
Length: 5cm, 7cm, 8cm, 9cm, 12cm.



BEVEL 45

Product Code:AK-LP-220

Available Diameter & Length
Diameter: 16ga, 15ga, 14ga.



BEVEL 30

Product Code:AK-LP-225

Available Diameter & Length
Diameter: 16ga, 15ga, 14ga.

Infusion Handles & (Cont.)



INFUSION HANDLE WITH ON/OFF CONTROL Model Number: **AK-LP-500**

Product Description:

A flow on/off control infusion handle for use with roller pump devices.

Handle diameter: **15/16" (Stainless Steel control)**

Valve: **Stainless Steel**

Female Luer Lock: **Stainless Steel**

Male Luer Lock: **Stainless Steel**



INFUSION HANDLE WITH ON/OFF CONTROL Model Number: **AK-LP-550**

Product Description:

A flow on/off control infusion handle for use with roller pump devices.

Handle diameter: **15/16" (control) 3/4" (handle)**

Valve: **Stainless Steel**

Female Luer Lock: **Stainless Steel**

Male Luer Lock: **Stainless Steel**



INFUSION HANDEL WITH ON/OF STOPCOCK CONTROL Model Number: **AK-LP-600**

Product Description:

Infusion Handle with Stopcock Control

Handle diameter: **3/4"**

Valve: **Brass with Electroless Nickel Plating**

Female Luer Lock: **Stainless Steel**

Male Luer Lock: **Stainless Steel**



INFUSION HANDEL WITH ON/OF STOPCOCK CONTROL Model Number: **AK-LP-650**

Product Description:

Infusion Handle with Stopcock Control

Handle diameter: **3/4"**

Valve: **Brass with Electroless Nickel Plating**

Female Luer Lock: **Stainless Steel**

Male Luer Lock: **Stainless Steel**



INFUSION HANDLE WITH SWIVEL LUER NUT

Model Number:

AK-LP-700

AK-LP-750 (for 1" Dia.)

Product Description:

Standard Infusion Handle, No Control with Swivel Nut

Handle diameter: **3/4"**

Female Luer Lock: **Stainless Steel**

Male Swivel Luer Lock: **Stainless Steel**



LUER LOCK ADJUSTABLE HANDLE

Model Number:

AK-LP-800(for small)

AK-LP-850(for large)

Product Description:

Luer Lock Adjustable Handle. Orientation between Luer Lock and Handle.

Handle diameter: **3/4" (for small) 1.0" (for large)**

Female Luer Lock: **Stainless Steel**

Male Luer Lock: **Stainless Steel**



Liposuction Cannula Catalogue

Infusion Handles & (Cont.)



INFUSION HANDLE WITH TUBE CONNECTOR

Model Number:
AK-LP-900
AK-LP-950(for 1" Dia.)

Product Description:

Infusion handle with tube connector for use with reusable infusion tubing.

Handle diameter: **3/4"**

Female Luer Lock: **Stainless Steel**

Barb: **Stainless Steel**



LUER LOCK ADJUSTABLE ASPIRATION HANDLE

Model Number: AK-LP-1000

Product Description:

Luer Lock Adjustable Handle ideal for rotating the

Luer Lock flat

Handle diameter: **3/4"**

Male Luer Lock: **Stainless Steel**



ASPIRATION HANDLE

Model Number:
AK-LP-1050(for large)
AK-LP-1100(for small)

Product Description:

Standard Aspiration Handle

Handle diameter: **3/4" (for small) 1.0" (for large)**

Male Luer Lock: **Stainless Steel**



INFUSION HANDLE WITH SWIVEL LUER NUT

Model Number:
AK-LP-1150
AK-LP-1200 L (for 1" Dia.)

Product Description:

Standard Infusion Handle, No Control with Swivel Nut

Handle diameter: **3/4"**

Female Luer Lock: **Stainless Steel**

Male Swivel Luer Lock: **Stainless Steel**



QUICK DISCONNECT CANNULA SET

Model Number:AK-LP-1250

Product Description:

A quick Disconnect Handle is a quick and easy interchangeable cannula hub.

Overall length: **6.65"**

Body: **All Stainless Steel**



QUICK LOCK CANNULASET

Model Number:AK-LP-1300

Product Description:

A quick Lock Handle is a quick and easy interchangeable Cannulas hub. Standard vented holes and swivel barb

Overall length: **6.65"**

Body: **Aluminum with type 2 blue anodized**

Lock : **Stainless Steel**



QUICK LOCK CANNULA SET

Model Number:AK-LP-1350

Product Description:

A quick Lock Handle is a quick and easy interchangeable Cannulas hub. Standard vented holes and swivel barb

Overall length: **6.65"**

Body: **Aluminum with type 2 blue anodized**

Lock : **Stainless Steel**



TOOMEY ASPIRATION CANNULA SET

Model Number:AK-LP-1400

Product Description:

Toomey Aspiration Handle

Overall length: **7.1"**

Body: **Delrin (Autoclavable)**



Liposuction Cannula Catalogue

Two Piece Cannuals & (Cont.)



LUER LOCK ADJUSTABLE ASPIRATION HANDLE

Model Number:AK-LP-1450

Product Description:

Luer Lock Adjustable Handle ideal for rotating the Luer Lock flat to the thumb rest

Handle diameter: **3/4"**

Overall length: **7.2"**

Body: **Aluminum with type 2 blue anodized**

Male Luer Lock: **Stainless Steel**



ASPIRATION HANDLE

Part Number:

AK-LP-1500 (for large)

AK-LP-1550 (for small)

Product Description:

Standard Aspiration Handle

Handle diameter: **3/4" (for small)**

1.0" (for large)

Overall length: **6.85"**

Body: **Aluminum with type 2 blue anodized**

Male Luer Lock: **Stainless Steel**

Injection Gun Set



INJECTION GUN SET 1ML, 3ML, 5ML, 10ML BD

Model Number:AK-LP-2000

Product Description:

Standard method of micro tissue injections.

All stainless steel

For 1ml, 3ml, 5ml, 10ml BD syringes

.15 ml per trigger on 1ml BD syringe

.32 ml per trigger on 3ml BD syringe

.62 ml per trigger on 5ml BD syringe

.88 ml per trigger on 10ml BD syringe



INJECTION GUN SET 3CC, 5CC, 10CC BD

Model Number:AK-LP-3000

Product Description:

Standard method of micro tissue injection. All stainless steel

For 3cc, 5cc, 10cc BD syringes

.32 ml per trigger on 3ml BD syringe

.62 ml per trigger on 5ml BD syringe

.88 ml per trigger on 10ml BD syringe



INJECTION GUN SET 3CC

Model Number:AK-LP-4000

Product Description:

Standard method of micro tissue injection. All stainless steel

For 3cc BD syringes

.32 ml per trigger on 3ml BD syringe



INJECTION GUN SET 5CC

Model Number:AK-LP-5000

Product Description:

Standard method of micro tissue injection.

All stainless steel

For 5cc BD syringes

.62 ml per trigger on 5ml BD syringe



Liposuction Cannula Catalogue

Injection Gun Set



INJECTION GUN SET 10CC

Model Number: AK-LP-6000

Product Description:

Standard method of micro tissue injection.
All stainless steel
For 10cc BD syringes
.88 ml per trigger on 10ml BD syringe



INJECTION GUN SET 20CC

Model Number: AK-LP-7000

Product Description:

Standard method of micro tissue injection.
All stainless steel
For 20cc BD syringes
1.5 ml per trigger on 20ml BD syringe



INJECTION GUN SET 30CC

Model Number: AK-LP-8000

Product Description:

Standard method of micro tissue injections. All stainless steel
For 30cc BD syringes
1.7 ml per trigger on 30ml BD syringe



INJECTION GUN SET 60CC

Model Number: AK-LP-9000

Product Description:

Standard method of micro tissue injections.
All stainless steel
For 60cc Monoject syringes
1.9 ml per trigger on 60ml Monoject syringe



Additional Syringe Accessories

| Part Number: | Product Description: |
|--------------|---------------------------|
| AK-LP-1500 | Plunger Adaptor Removal |
| AK-LP-1550 | Plunger adaptor for 3 ml |
| AK-LP-1560 | Plunger adaptor for 5 ml |
| AK-LP-1570 | Plunger adaptor for 10 ml |
| AK-LP-1580 | Plunger adaptor for 20 ml |
| AK-LP-1590 | Plunger adaptor for 30 ml |
| AK-LP-1600 | Plunger adaptor for 35 ml |
| AK-LP-1650 | Plunger adaptor for 60 ml |

Product Description:

This additional accessories are ideal for a surgeon who wants to have two or more syringes ready for tissue re-injection procedure.



Liposuction Cannula Catalogue

Toomey & Accessories



TOOMEY TO LUER LOCK SYRINGE

Model Number:AK-LP-2002

Product Description:
Toomey 60cc Monoject to Luer Lock Syringe
Stainless Steel and Stainless Steel hub



TOOMEY TO LUER LOCK SYRINGE

Model Number:AK-LP-2003

Product Description:
Toomey 60cc Monoject to Luer Lock Syringe
Female Luer Lock made of Stainless Steel and Aluminum
hub



TOOMEY TO LUER LOCK NEEDLE

Model Number:AK-LP-2004

Product Description:
Toomey 60cc Monoject to Luer Lock Needle
Stainless Steel hub



TOOMEY TO LUER LOCK NEEDLE

Model Number:AK-LP-1450

Product Description:
Toomey 60cc Monoject to Luer Lock Needle
Male Luer Lock made of Stainless Steel and Aluminum
hub



SYRINGE BARB ADAPTOR

Part Number:
AK-LP-1460 SBA60SS
AK-LP-1465 SBA35SS
AK-LP-1462 SBA10SS

Product Description:
Syringe Barb Adapter
Made of Stainless Steel



LUER LOCK TO LUER LOCK

Model Number:AK-LP-1468

Product Description:
Luer Lock syringe to Luer Lock syringe soft Tissues
Transfer system
Made of Stainless Steel



Sterilization Tray



STERILIZATION TRAY

Model No AK-LP-3600

Trays are sterilizable by all standard methods of sterilization. Trays are designed to encase instruments during processing. Safe for all standard sterilization methods. Autoclave/dry heat up to 360°F, ETO, and Cold Solutions.

Flexibility of layout. Finger tip mats.

Unique, "grid" system in all bases, trays, and lids makes it easy to install mats or inserts to protect your instruments.

4.0x7.5x1.5 Micro Tray Base Lid & Mat



Commissioner Road, Street # 4.
Muhammad Pura.
Sialkot-51310. Pakistan

Phone: +92 332 04444 57
E-Mail: info@aksurgicalpro.com
Website: www.aksurgicalpro.com